

SIEMENS MAGNETOM Prisma syngo MR D13D

\\USER\NIMH\Hebart\localizer\AAHead_Scout
 TA:0:17 PAT:2 Voxel size:1.6×1.6×1.6 mm Rel. SNR:1.00 :fl

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A10.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
TE	1.37 ms
Averages	1
Concatenations	1
Filter	Normalize
Coil elements	HC1-7
AutoAlign	Head

Contrast

Flip angle	8 deg
Averaging mode	Short term
Measurements	1
Reconstruction	Magnitude

Resolution

Base resolution	160
Phase resolution	100 %
Phase partial Fourier	6/8
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Prescan Normalize	Off
Normalize	On
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	69 %
Slice partial Fourier	6/8

Geometry

Nr. of slab groups	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A10.0 H0.0 mm
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	128
Multi-slice mode	Sequential
Series	Ascending
Nr. of sat. regions	0
Position mode	L-P-H
Special sat.	None
Table position	P
Inline Composing	Off

System

Body	Off
HC1	On
HC3	On
HC5	On
NC1	Off
HC2	On
HC4	On
HC6	On
HC7	On
NC2	Off
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
Coil Select Mode	Off - AutoCoilSelect
B0 Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.255397 MHz
Correction factor	1
SRFExcit 1H	39.852 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	3D
Averaging mode	Short term
Multi-slice mode	Sequential
Asymmetric echo	Weak
Contrasts	1
Bandwidth	540 Hz/Px
RF pulse type	Fast
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HC1-7
Acquisition duration	0 ms
Mode	Off

BOLD

Time to center	7.5 s
Subtract	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
Contrasts	1
Save original images	On

SIEMENS MAGNETOM Prisma syngo MR D13D

\\USER\NIMH\Hebart\localizer\T1w_MPR

TA:6:38 PAT:2 Voxel size:0.8x0.8x0.8 mm Rel. SNR:1.00 :tfl

Properties

Prio Recon	Off
Load to viewer	Off
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	23.1 %
FoV read	256 mm
FoV phase	93.8 %
Slice thickness	0.80 mm
TR	2400.0 ms
TE	2.24 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

Contrast

Magn. preparation	Non-sel. IR
TI	1060 ms
Flip angle	8 deg
Fat suppr.	Water excit. fast
Water suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement

Resolution

Base resolution	320
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	On
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	Off

Geometry

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	23.1 %
Slices per slab	208
Multi-slice mode	Single shot
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Water excit. fast
Water suppr.	None
Special sat.	None
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Coil Select Mode	Off - All
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
F >> H	256 mm
A >> P	240 mm
R >> L	167 mm
Frequency 1H	123.255397 MHz
Correction factor	1
ExcitWEs 0 1H	59.778 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	Non-sel. IR
TI	1060 ms
Dark blood	Off
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Long term
Multi-slice mode	Single shot
Reordering	Linear
Asymmetric echo	Allowed
Bandwidth	210 Hz/Px
Flow comp.	No
Echo spacing	8.1 ms
Turbo factor	256
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Subtract	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
Save original images	On

SIEMENS MAGNETOM Prisma syngo MR D13D

\\USER\NIMH\Hebart\localizer\pRF
TA:7:51 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	60
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	P >> A
AutoAlign	---
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	1500 ms
TE	33.00 ms
Multi-band accel. factor	3
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	75 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	308
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	96
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	60
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	P >> A
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Coil Select Mode	Off - All
B0 Shim mode	Advanced
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	180.00 deg
R >> L	192 mm
A >> P	192 mm
F >> H	120 mm
Frequency 1H	123.255397 MHz
Correction factor	1
MBExc 1H	366.626 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Contrasts	1
Bandwidth	2264 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	0.55 ms
EPI factor	96
Gradient mode	Performance
Excitation	Standard
RF spoiling	Off
Online multi-band recon.	Online
Physio recording	Off
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	Off
Temp. highpass filter	Off
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off
Contrasts	1

SIEMENS MAGNETOM Prisma syngo MR D13D

\\USER\NIMH\Hebart\localizer\6cat_localizer
TA:4:39 PAT:Off Voxel size:2.0×2.0×2.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	60
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	P >> A
AutoAlign	---
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	1500 ms
TE	33.00 ms
Multi-band accel. factor	3
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	75 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	180
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	96
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	60
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	P >> A
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Coil Select Mode	Off - All
B0 Shim mode	Advanced
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	180.00 deg
R >> L	192 mm
A >> P	192 mm
F >> H	120 mm
Frequency 1H	123.255397 MHz
Correction factor	1
MBExc 1H	366.626 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Contrasts	1
Bandwidth	2264 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	0.55 ms
EPI factor	96
Gradient mode	Performance
Excitation	Standard
RF spoiling	Off
Online multi-band recon.	Online
Physio recording	Off
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	Off
Temp. highpass filter	Off
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off
Contrasts	1

SIEMENS MAGNETOM Prisma syngo MR D13D

\\USER\NIMH\Hebart\localizer\T2w_SPC
TA:5:57 PAT:2 Voxel size:0.8x0.8x0.8 mm Rel. SNR:1.00 :spc

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	0.0 %
FoV read	256 mm
FoV phase	93.8 %
Slice thickness	0.80 mm
TR	3200 ms
TE	564.0 ms
Concatenations	1
Filter	Prescan Normalize, Image Filter
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Fat suppr.	None
Restore magn.	Off
Measurements	1
Reconstruction	Magnitude
Multiple series	Each measurement

Resolution

Base resolution	320
Phase resolution	100 %
Phase partial Fourier	Allowed
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Integrated
Image Filter	On
Intensity	Sharp
Edge Enhancement	3
Smoothing	3
Unfiltered images	Off
Distortion Corr.	Off
Accel. factor 3D	1
Unfiltered images	On
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	Off

Geometry

Nr. of slab groups	1
Slabs	1
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	208
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Special sat.	None
Special sat.	None
Table position	P
Inline Composing	Off
Restore magn.	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	Head > Brain
Coil Select Mode	Off - All
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
F >> H	256 mm
A >> P	240 mm
R >> L	167 mm
Frequency 1H	123.255397 MHz
Correction factor	1
SPC_nsExc 1H	224.167 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Trigger delay	0 ms
Magn. preparation	None
Dark blood	Off
Resp. control	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Bandwidth	744 Hz/Px
Flow comp.	No
Allowed delay	0 s
Echo spacing	3.86 ms
Adiabatic-mode	Off
Turbo factor	314
Echo train duration	1166 ms
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
Flip angle mode	T2 var
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Organ under exam.	Standard
Tissue T1	940 ms
Tissue T2	100 ms

BOLD

Subtract	Off
StdDev	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
Save original images	On

SIEMENS MAGNETOM Prisma syngo MR D13D

\\USER\NIMH\Hebart\localizer\AP reference
TA:6:09 PAT:Off Voxel size:2.0x2.0x2.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	60
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	1500 ms
TE	33.00 ms
Multi-band accel. factor	3
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Magn. preparation	None
Flip angle	75 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	240
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	96
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	60
Dist. factor	0 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Coil Select Mode	Off - All
B0 Shim mode	Advanced
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	192 mm
A >> P	192 mm
F >> H	120 mm
Frequency 1H	123.255397 MHz
Correction factor	1
MBExc 1H	366.626 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Magn. preparation	None

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	Off
Averaging mode	Long term
Multi-slice mode	Interleaved
Contrasts	1
Bandwidth	2264 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	0.55 ms
EPI factor	96
Gradient mode	Performance
Excitation	Standard
RF spoiling	Off
Online multi-band recon.	Online
Physio recording	Off
Triggering scheme	Standard
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	Off
Temp. highpass filter	Off
Threshold	4.00
Paradigm size	3
Motion correction	Off
Spatial filter	Off
Delay in TR	0 ms
Distortion Corr.	Off
Contrasts	1

SIEMENS MAGNETOM Prisma syngo MR D13D

\\USER\NIMH\Hebart\localizer\T2STAR SAG 3DEPI
TA:5:48 Voxel size:0.7×0.7×0.7 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	R6.9 P7.9 F5.5 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Phase oversampling	0 %
Slice oversampling	0.0 %
FoV read	250 mm
FoV phase	81.3 %
Slice thickness	0.65 mm
TR	64 ms
TE	35.0 ms
Averages	1
Concatenations	1
Filter	Raw filter, Prescan Normalize
Coil elements	HEA;HEP

Contrast

MTC	Off
Flip angle	10 deg
Fat suppr.	Water excit. normal
Averaging mode	Long term
Measurements	1
Reconstruction	Magn./Phase
Multiple series	Each measurement

Resolution

Base resolution	384
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Distortion Corr.	Off
Hamming	Off
Unfiltered images	Off
Prescan Normalize	On
Raw filter	On
Intensity	Strong
Slope	64
Elliptical filter	Off
Slice resolution	100 %
Slice partial Fourier	Off

Geometry

Nr. of slab groups	1
Slabs	1
Dist. factor	50 %
Position	R6.9 P7.9 F5.5 mm
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	256
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Water excit. normal
Special sat.	None
Special sat.	None
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	Head > Basis
Coil Select Mode	Default
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	R6.9 P7.9 F5.5 mm
Rotation	0.00 deg
F >> H	250 mm
A >> P	204 mm
R >> L	167 mm
Frequency 1H	123.255397 MHz
Correction factor	1
EPI_WEWEs 0 1H	48.905 V
Gain	Low
Table position	0 mm
Img. Scale. Cor.	1.000

Physio	1st Signal/Mode	None
Inline	Inline Composing	Off
	Distortion correction	Off
Sequence	Introduction	On
	Dimension	3D
	Averaging mode	Long term
	Multi-slice mode	Interleaved
	Bandwidth	394 Hz/Px
	Free echo spacing	Off
	Echo spacing	2.78 ms
	EPI factor	15
	RF pulse type	Normal
	Gradient mode	Fast
	RF spoiling	On
	TX/RX delta frequency	0 Hz
	TX Nucleus	None
	TX delta frequency	0 Hz
	Coil elements	HEA;HEP
	Acquisition duration	0 ms
BOLD	Subtract	Off
	StdDev	Off
	MIP-Sag	Off
	MIP-Cor	Off
	MIP-Tra	Off
	MIP-Time	Off
	Save original images	On
	Distortion Corr.	Off
	Save original images	On

SIEMENS MAGNETOM Prisma syngo MR D13D

\\USER\NIMH\Hebart\localizer\TOF_3D_multi-slab TA:8:18 PAT:2 Voxel size:0.3×0.3×0.5 mm Rel. SNR:1.00 :fl_r

Properties	Prio Recon	Off
	Load to viewer	On
	Inline movie	Off
	Auto store images	On
	Load to stamp segments	On
	Load images to graphic segments	Off
	Auto open inline display	Off
	Wait for user to start	Off
	Start measurements	single

Routine

Nr. of slab groups	1
Slabs	6
Dist. factor	-20.00 %
Position	L0.0 P10.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	20.0 %
FoV read	200 mm
FoV phase	90.6 %
Slice thickness	0.50 mm
TR	21.0 ms
TE	3.43 ms
Averages	1
Concatenations	6
Filter	Normalize
Coil elements	HEA;HEP

Contrast

TD	0.000 ms
MTC	Off
Flip angle	18 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Short term
Measurements	1
Reconstruction	Magnitude

Resolution

Base resolution	384
Phase resolution	95 %
Phase partial Fourier	Off
Interpolation	On
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
TD	0.000 ms
Accel. factor 3D	1
Prescan Normalize	Off
Normalize	On
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
POCS	Off
Slice resolution	50 %
Slice partial Fourier	7/8

Geometry

Nr. of slab groups	1
Slabs	6
Dist. factor	-20.00 %
Position	L0.0 P10.0 H0.0 mm
Phase enc. dir.	R >> L
Phase oversampling	0 %
Slice oversampling	20.0 %
Slices per slab	40
Multi-slice mode	Sequential
Series	Ascending
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Water suppr.	None
Special sat.	Tracking F
Gap	10 mm
Thickness	40 mm
Special sat.	Tracking F
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect
B0 Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm
Frequency 1H	123.255397 MHz
Correction factor	1
SRFExcit 1H	295.639 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

1st Signal/Mode	None
Dark blood	Off

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Averaging mode	Short term
Multi-slice mode	Sequential
Asymmetric echo	Allowed
Contrasts	1
Bandwidth	186 Hz/Px
Flow comp.	Yes
Gradient mode	Fast
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

TONE ramp	70 %
Flow direction	H >> F
3D centric reordering	On
Subtract	Off
StdDev	Off
MIP-Sag	On
MIP-Cor	On
MIP-Tra	On
MIP-Time	Off
Save original images	On
Distortion Corr.	Off
Contrasts	1
Save original images	On

SIEMENS MAGNETOM Prisma syngo MR D13D

\\USER\NIMH\Hebart\localizer\gre_field_mapping_2mm
TA:1:09 Voxel size:3.0×3.0×3.0 mm Rel. SNR:1.00 :fm_r

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	49
Dist. factor	0 %
Position	L0.3 P13.8 F33.9 mm
Orientation	Transversal
Phase enc. dir.	P >> A
AutoAlign	---
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	520.0 ms
TE 1	4.92 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	HEA;HEP

Contrast

MTC	Off
Flip angle	60 deg
Fat suppr.	None
Averaging mode	Long term
Measurements	1
Reconstruction	Magn./Phase
Multiple series	Off

Resolution

Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	49
Dist. factor	0 %
Position	L0.3 P13.8 F33.9 mm
Phase enc. dir.	P >> A
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	None
Special sat.	None
Special sat.	None
Table position	P
Inline Composing	Off

System

Body	Off
HEP	On
HEA	On
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	L0.3 P13.8 F33.9 mm
Rotation	180.00 deg
R >> L	192 mm
A >> P	192 mm
F >> H	147 mm
Frequency 1H	123.255397 MHz
Correction factor	1
01GreFCE 1H	127.955 V
Gain	High
Table position	0 mm
Img. Scale. Cor.	1.000

Physio

Inline

Inline Composing	Off
Distortion correction	Off

Sequence

Introduction	On
Dimension	2D
Averaging mode	Long term
Multi-slice mode	Interleaved
Asymmetric echo	Off
Contrasts	2
Bandwidth	300 Hz/Px
Flow comp.	Yes
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HEA;HEP
Acquisition duration	0 ms
Mode	Off

BOLD

Distortion Corr.	Off
Contrasts	2